

Woodlyn® Millennium Series Stand



The Woodlyn Millennium Instrument Stand was designed to offer state-of-the-art functionality for the operator while conveniently accommodating patients. The control console is angled for easy access, and features two illuminated chair switches (up and down), overhead light illumination adjustment, and three rechargeable wells with indicators. There is a low voltage power supply (for optional BIO) which has a built in hanger. When the instrument is placed on the hanger, power is shut off. When in use, the power is automatically turned on. Illumination control for the BIO is conveniently positioned on the face of the control console. The slit lamp arm was designed for maximum travel; this allows taller patients to be examined with ease, and also allows the arm to go low enough to reach patients in a wheelchair, without reseating the patient.

SPECIFICATIONS

Stand

Weight of base unit and OHL: 400 lbs.
Column Diameter: 1 3/4"

Overhead Lamp

20W Halogen

Lower Arm

Height in lowest position at instrument end: 23 1/2"
Height in highest position at instrument end: 37"

Counter Balance

Factory counter balanced for an instrument weighing 38 lbs. Additional weights available to counter balance instruments up to 50 lbs.

Console Controls

Overhead lamp rheostat
Auxiliary corded instrument rheostat (BIO)
Chair Functions: Power up and down
Three rechargeable wells